


IMPLANT MADE OF BLACKARMOR®

A nonmetallic, radiolucent material
for reduced imaging artifacts

ANATOMICAL SHAPE

For precise and secure implantation



SECURE FIXATION TO BONE

Rough titanium coating

inspired by nature – built by icotec

ICOTEC CERVICAL CAGE

Factsheet




EXCELLENT IMAGING QUALITY

focused on bone – not metal

ICOTEC CERVICAL CAGE

icotec Cervical Cages – Anatomic

Description	Depth	Width	Height	Reference number
icotec Cervical Cage, anatomic	13 mm	13 mm	5 mm	13-41-0513
			6 mm	13-41-0613
			7 mm	13-41-0713
icotec Cervical Cage, anatomic	13 mm	16 mm	5 mm	13-41-0516
			6 mm	13-41-0616
			7 mm	13-41-0716
			8 mm	13-41-0816
icotec Cervical Cage, anatomic	13 mm	18 mm	5 mm	13-41-0518
			6 mm	13-41-0618
			7 mm	13-41-0718
			8 mm	13-41-0818



Information needed to use the device is made available in electronic form; current and previous versions can be downloaded in electronic form at ifu.icotec-medical.com (code = [REFI](#)) or can be requested by email or phone from icotec. On request, icotec will provide a paper version within seven calendar days at no charge. The electronic versions can be viewed with a freely available PDF reader (e.g., Adobe Acrobat Reader, which can be downloaded at www.adobe.com).



ifu.icotec-medical.com

ico | cæt

 icotec ag

Industriestrasse 12
 9450 Altstätten
 Switzerland
 Phone: +41 71 757 00 00
 Fax: +41 71 757 00 01
info@icotec.ch
www.icotec-medical.com

013005.2019-01