


BLACKARMOR® CARBON/PEEK RODS

Product Overview

 Swiss Quality Product



Carbon/PEEK implant material from icotec



| Description | icotec Carbon/PEEK rod straight | icotec Carbon/PEEK rod curved | icotec Carbon/PEEK rod curved | icotec Carbon/PEEK J-rod: combination of straight and curved |
|------------------|-------------------------------------|------------------------------------|-------------------------------------|--------------------------------------------------------------|
| Diameter | 5.5 mm | 5.5 mm | 5.5 mm | 5.5 mm |
| Radius | 000 | 120 | 450 | 000 & 450 |
| Length | 40 - 100 mm in 10 mm increments | 40 - 100 mm in 10 mm increments | 100 - 160 mm in 20 mm increments | 60/70 mm |
| | 100 - 160 mm in 20 mm increments | | | 50/100 mm |
| Reference Number | 16-55-000xx | 16-55-012xx | 16-55-045xx | 60/100 mm |
| | | | | 16-55-00607-9 |
| | | | | 16-55-00510-9 |
| | | | | 16-55-00610-9 |

Initial rod portfolio. Additional rods of intermediate lengths can be made available upon demand.