

Factsheet

icotec Cervical Cage

Anatomical shape

For precise and secure implantation




Implant made of BlackArmor®

A nonmetallic, radiolucent material for reduced imaging artifacts



Secure fixation to bone

Rough titanium coating



Excellent imaging quality


Focused on bone – not metal



icotec Cervical Cages – Anatomic

Description	Depth	Width	Height	Reference number
icotec Cervical Cage, anatomic	13 mm	13 mm	5 mm	12-41-0513* 13-41-0513
			6 mm	12-41-0613* 13-41-0613
			7 mm	12-41-0713* 13-41-0713
icotec Cervical Cage, anatomic	13 mm	16 mm	5 mm	12-41-0516* 13-41-0516
			6 mm	12-41-0616* 13-41-0616
			7 mm	12-41-0716* 13-41-0716
			8 mm	12-41-0816* 13-41-0816
icotec Cervical Cage, anatomic	13 mm	18 mm	5 mm	12-41-0518* 13-41-0518
			6 mm	12-41-0618* 13-41-0618
			7 mm	12-41-0718* 13-41-0718
			8 mm	12-41-0818* 13-41-0818

* New reference number as of 01/2024 (without any modification to the implant).



icotec ag

Contact in the USA:

222 Pitkin Street
East Hartford CT 06108

United States
Phone (860) 404-6999
info@icotec-medical.com
www.icotec-medical.com